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Muffel-System thermopress 400

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Instructions for use

thermopress 400 muffle system

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage.

1. Product description

The thermopress 400 muffle system (REF 110 0040 4) is a transfer tool made of stainless steel and consists of 5 components. It is used to be able to insert investment material muffles (rings) into thermopress 400 injection moulding units and to process them. The user is enabled to model, sprue and invest wax models in accordance with familiar dental procedures. The scope of delivery includes an additional conversion kit (pressure plate with threaded bar), which enables the user to process the muffle system with older thermopress 400 units.

The scope of delivery includes the following components:

- thermopress 400 metal frame (REF 360 0128 7)
- thermopress 400 base plate (REF 360 0128 6)
- thermopress 400 base former (REF 360 0128 5)
- thermopress 400 silicone sleeve with O-ring (REF 360 0128 0)
- thermopress 400 conversion kit for older units (REF 360 0128 8)

2. Indication range

The thermopress 400 muffle system is used for processing the thermoplastic materials Bio XS and bre.dentan HP. Thanks to individual expansion control of the investment material (Brevest Rapid 1) fast and simple fabrication of fixed dental restorations (crowns and bridges) with the biocompatible resins (Bio XS and bre.dentan HP) of the thermopress system is enabled.

3. Contraindications

The thermopress 400 muffle system is not suitable for processing other thermoplastic materials, such as bre.crystal HP, Bio Dentaplast, bre.flex and bre.flex 2nd edition. bredent is not able to provide any processing information on processing other thermoplastic materials from competitors either. Usage of the muffle system in other resin injection systems is at the user's own risk.

4. Safety recommendations and hazard warning

Note:

The muffle needs to be preheated to 150°C for 45 minutes before the injection process. As a result, there is a risk of getting burned when taking the muffle out of the preheating furnace and placing it into the thermopress 400 unit.
Please use only suitable protective gloves or muffle (furnace) tongs.
Please observe also the relevant safety recommendations when processing phosphate-bonded investment materials (see instructions for use for Brevest Rapid 1).

5. Durability and storage

No special recommendations need to be observed!

6. Processing/Use

6.1 Modellring

The wax model is fabricated in the usual way. A minimum wall thickness of 0.5 mm must be ensured. Chamber preparation by the dentist is required (see fig. 1).

6.2 Spruing

Parts/components required:

- Base former, silicone sleeve
- Beels of wax pattern, 2.0 mm (REF 430 0116 0)
- Reels of wax pattern, 5.0 mm (REF 430 0118 5)
- thermopress acrylic sprue wax, 10 mm (REF 430 0741 0)

The resin injection channel (sprue – 10 mm) is placed in the opening of the base former and cut off to obtain a length of 10 mm. A wax channel (sprue – 2 mm) is attached to each bridge unit (vertical to the occlusion surface). This 2 mm wax sprue is affixed directly to the incisal edge of anterior teeth. A 5 mm injection bar is fastened horizontally to each 2 mm wax sprue at a distance of 3 mm. One end is attached to the main sprue and the other end is bent upwards to serve as a pressure compensation sprue. A hot wax knife is used to reduce the diameter of this pressure compensation sprue (half of its original size) (see fig. 2). The excess material can escape during the injection process and the object will not be damaged by excessive pressure. The bottom plate features a recess to take up the excess material.

Important:

Sufficient distance (at least 5 mm) between the wax model and the pressure compensation sprues towards the inner wall of the silicone sleeve must be ensured (see fig. 3).

6.3 Preparation of the muffle

Parts/components/devices required:


- Base former, metal frame, base plate
- Silicone sleeve
- Isosil (REF 520 IS12 5)
- Brevest Rapid 1 (REF 570 160R 4)
- ecovac vacuum mixing unit (REF 140 0093 0)
- Dental vibrator
- Preheating furnace
- Expando-Rock (REF 570 OERS 5)

Once the wax model has been attached to the base former, Isosil is sprayed onto the inner side of the silicone sleeve and the base former. Then the silicone sleeve is placed on the base former. Make sure to keep a minimum distance of 5 mm between the wax model and the silicone sleeve (see fig. 3). 160 g of Brevest Rapid investment material need to be mixed with a 82 % concentration of Brevest R in a vacuum mixing unit for 60 sec and the formation of bubbles must be avoided. Then the investment material is filled into the muffle (use a dental vibrator). After a hardening time of 20 min., the muffle (without the silicone sleeve and the base former) is preheated in the preheating furnace for 60 min. at 900°C and the wax is eliminated. Then the muffle needs to cool down to room temperature.

Afterwards the base plate is placed into the metal frame and this way the bottom is sealed tightly (see fig. 4).

The investment material muffle that has cooled down to room temperature is centered on the base plate in this metal frame (see fig. 5). The base former also serves as a pressure plate and is placed on the investment material muffle using the guide pin.

The empty space between the investment material muffle and the metal frame is filled with Expando-Rock (26 ml concentrate : 100 g powder) up to the upper edge of the metal frame (see fig. 7).

Please make sure to cover the injection opening for a short time.

6.4 Injection process

After filling the investment material muffle in the metal frame, the plaster needs to harden approx. 30 min. Then the entire muffle is preheated at 150°C for 45 min. At the same time, the program of the injection moulding system thermopress 400 for processing Bio XS or bre.dentan HP can already be selected and started. The heating process starts. Once the material to be used has melted, which is indicated by an acoustic signal, the metal muffle is removed from the furnace and placed in the thermopress 400 unit; first the front clamping screw (without clamping effect) and then the lateral clamping screw are hand-tightened. Then both clamping screws are tightened further. Afterwards the injection process can be started.

Important:

The front clamping device of older units features a thicker pressure plate (see fig. 8). This pressure plate needs to be exchanged using the conversion kit enclosed to the muffle set. New units already have a thin pressure plate which does not need to be exchanged.

The unit is operated in accordance with the operating instructions.

6.5 Devesting

Once the press plunger of the thermopress 400 unit has returned to the initial position, both clamping screws can be loosened and the lateral bracket needs to be lifted up to avoid damaging the unit when ejecting the cartridge.

After ejecting the cartridge, the muffle can be removed from the thermopress 400 unit. Use a plastic hammer to tap lightly onto the base plate at the bottom of the metal frame to remove the entire investment material muffle with the surrounding Expando-Rock material from the metal frame. The injected object is carefully removed using a pneumatic deflasking chisel or light hammer taps. Investment material residues are removed from the injected object using 125 µm Al₂O₃ and a pressure of 3 bar. The injection sprues can be removed with a diamond disc.

7. Technical data

Weight (metal frame) = 2 kg

Total weight of muffle system = 3 kg

Dimensions of the metal frame – w x h x d in mm = 67 x 127 x 84

8. Symbols

 Observe instructions for use

9. Additional information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend to read the instructions for use again before using a new package. These non-binding user recommendations are based on our own experience. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since bredent does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.
The product may only be used by dentists, dental technicians and adequately trained dental staff. Only original tools and components may be used for processing.

If you have questions concerning handling, please call **+49 – 73 09 / 8 72-22!**

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